

Surgical Tissue Management System
Ultrasonic Generator

**USG-400** 



### Main features

# USG-400 supports

- THUNDERBEAT, Exclusively delivering both ultrasonic and bipolar energy for reliable vessel sealing and the fastest in class tissue cutting, when combined with ESG-400, electrosurgical generator.
- SONICBEAT, Delivering ultrasonic energy for the simultaneous cut and coagulation of tissue.

## USG-400 Features and Benefits

- Graphical User Interface for intuitive operation, ease of use, and rapid troubleshooting.
- Instrument recognition and automatic application of default settings for ease of use.
- Hand switch and foot switch activation options to meet each surgeon's individual preferences.
- Automatic mist and smoke evacuation to maintain a clear laparoscopic view reducing delays associated with compromised visualization, when combined with UHI series insufflator.
- A dedicated cart for transportation and storage.





USG-400 and ESG-400 (optional) installed on TC-E400 (optional)

#### **USG-400 Specifications**

Power Supply	Voltage	120/220 - 240 V AC
	Voltage fluctuation	Within ±10%
	Frequency	50/60 Hz
	Frequency fluctuation	Within ±1 Hz
	Input	360 VA
	Fuse rating	4 A, 250 V
	Fuse size	ø 5 × 20 mm
Size	Dimensions	375 (W) × 156 (H) × 480 (D) mm
	Weight	9 kg
Classification (medical electrical equipment)	Type of protection against electric shock	Class I
	Degree of protection against electric shock	CF type
	Degree of protection against explosion	The ultrasonic generator should be kept away from flammable gases.
Output	Applicable instrument	THUNDERBEAT*1
		SONICBEAT
	ON/OFF operation	Footswitch
		Handswitch
	Output setting	Setting possible in the Set Screen displayed on the touch panel.
	Frequency	47.0 kHz

<sup>\*1</sup> When combined with ESG-400, electrosurgical generator.

## TC-E400 Energy Cart Specifications

Dimensions	515 (W) × 1047 (H) × 715 (D) mm
Shelf Width	378 mm
Weight	38 kg
Load Capacity	Equipment shelf 40 kg*2
(when equally distributed)	Foot switch shelf 20 kg
(when equally distributed)	Drawer 5 kg
Casters	4 x ø 125 mm dia., with brakes on all the casters.

<sup>\*2</sup> Do not place an instrument with a larger footprint than the shelf plate on the EPF Cart.

### Ordering information

USG-400	Ultrasonic Generator	
Option		
ESG-400	Electrosurgical Generator *3	
TC-E400	Energy cart for USG-400 and ESG-400	
MAJ-1869	Footswitch for SONICBEAT	
MAJ-1870	Footswitch for THUNDERBEAT	
MAJ-1871	Communication cable short *3	
MAJ-1872	Communication cable long *4	
MAJ-1873	Adapter for the connection to UHI-2 or UHI-3 insufflator *4	
MAJ-1876	Docking fixture *3	

<sup>\*3</sup> Required for the use of THUNDERBEAT

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer



<sup>\*4</sup> Required for the automatic mist and smoke evacuation function